

## TECHNICAL DATA SHEET



### Sterile 70% Isopropanol

- Gamma-irradiated to a sterility level of  $10^6$  with independent QC audits to ensure sterility
- Each lot tested for endotoxins
- Lot traceable
- Lot-specific information on each shipment simplifies record keeping:
  - › Certificate of Processing confirming radiation dosage
  - › Certificates of Compliance and Analysis
  - › Expiration date marked on each package

### Description

Texwipe's Sterile 70% Isopropanol contains 70% by volume USP-grade Isopropanol (isopropyl alcohol) and 30% USP-purified water. This solution is submicron filtered, filled into cleaned containers, double-bagged and gamma-irradiated. Each bottle is lot coded with an expiration date for easy recordkeeping.

Bottles are fully-assembled and ready for use – no unpacking, handling or attaching spray heads is required. The adjustable spray head offers stream delivery or coarse spray.

### Applications

- ▶ Surface cleaning
- ▶ Use on gloved hands in sterile suites
- ▶ Wipe down for pass through to controlled environments
- ▶ Ideally suited for use with sterile cleanroom wipers

### Features

- ▶ 0.2  $\mu$ m filtered 70% USP-grade Isopropanol (Isopropyl Alcohol)
- ▶ Non-atomizing adjustable spray
- ▶ 5 color indicator dots (for identification)
- ▶ Double-bagged in solvent-safe packaging
- ▶ Meets FDA guideline for medical device bacterial endotoxin limit of not more than 20 EU/device

### Benefits

- ▶ Ready-to-use sprayer eliminates handling and assembly, maintaining sterility
- ▶ Removes many surface contaminants
- ▶ Evaporates and leaves extremely low residue
- ▶ Prefiltered and sterilized for pass through and use in critical environments
- ▶ Eliminates mixing, filtration, processing and QC expense of in-house production of sterile IPA

### Products

Number	Description	Size	Packaging
<i>Mini-Cases</i>			
<b>TX3274</b>	Sterile 70% Isopropanol trigger spray	16 fluid ounces (473 mL)	4 polybottles/case
<i>Full-Cases</i>			
<b>TX3270</b>	Sterile 70% Isopropanol trigger spray	16 fluid ounces (473 mL)	12 polybottles/case
<b>TX3273</b>	Sterile 70% Isopropanol trigger spray	32 fluid ounces (946 mL)	12 polybottles/case
<b>TX8270</b>	Sterile 70% Isopropanol trigger spray	8 fluid ounces (237 mL)	12 polybottles/case

**North America:** 1210 South Park Drive, Kernersville, NC 27284 • Tel 336.996.7046 • Fax 336.996.2297  
**Europe/Middle East:** Saffierlaan 5, 2132 VZ Hoofddorp, The Netherlands • Tel +31 88 1307 400 • Fax +31 88 1307 499  
**Asia/Pacific:** ITW Singapore (Pte) Ltd., 500 Chai Chee Lane, Singapore 469024 • Tel +65 6468 9433 • Fax +65 6468 6772

---

### Product Characteristics

---

Property	Description
IPA Quality	Meets the requirements of USP Isopropyl Alcohol
Water Quality	Meets the requirements of USP Purified Water
Appearance	Clear, colorless liquid
Odor	Characteristic alcohol
% IPA Content	68% – 72%
Specific Gravity	0.872 – 0.883
Acidity	Not more than 1.0 mL of 0.020 N Sodium Hydroxide
Limit of NVR	Does not exceed 5 mg (0.01%)

---

**Certificate of Analysis provided with every lot.**

**Texwipe Regulatory Certificates**

Texwipe is the only wiper company to be registered: ISO 9001 • ISO 13485 • ISO 14001 • OHSAS 18001