

50  
YEARS




 **Texwipe®**  
An ITW Company

# Cleanroom Swabs

Microdenier | Polyester | Foam | ESD | Cotton | Specialty

Cleaning Validation & TOC Validation Kits | Alpha® | Absorbond® | CleanFoam®  
ESD-Safe | Cotton | Specialty | General Purpose

# Table of Contents




**Cleaning Validation**  
**Pg. 1**



**Cleaning Validation Kits**  
**Pg. 1**



**Microdenier**  
**Pg. 2**



**Polyester**  
**Pg. 3**




**CleanFoam®**  
**Pg. 4**



**ESD-Safe**  
**Pg. 6**



**Cotton**  
**Pg. 7**



**Specialty**  
**Pg. 7**



**General Purpose**  
**Pg. 8**



**Solvent Compatibility  
Guide / INDEX**  
**Pg. 9**

Our highly professional team will help you select and develop products for any critical environment application. Contact information is located on the back of this brochure or visit us at [www.texwipe.com](http://www.texwipe.com).

## About Texwipe

Texwipe's commitment to innovation, leadership and quality in cleanroom consumable products spans more than fifty years. We invest in technology to respond to our customers' evolving needs in contamination control.

### *Texwipe's Core Values are:*

- **Innovation** - Texwipe pioneers the latest technologies to provide innovation in contamination control products and processes.
- **Quality** - Texwipe product quality is maintained by the most advanced testing and quality control standards in the industry.
- **Technology Leadership** - Texwipe leads our industry in testing metrology, methods and processes to reduce contamination.

Throughout Texwipe's global operations, we support our customers with products designed to exceed the requirements for cleanroom consumable products. Our highly professional team will help you select and develop products for any critical environment application.



## Texwipe CleanTips®

Texwipe Cleanroom Swabs, with their trademark green handles, are the standard in cleanrooms across the globe. Engineered for superior cleanliness and consistency for a variety of applications, Texwipe has a swab that will provide the best contamination performance available.

Our products range from Microdenier, Polyester, CleanFoam®, ESD-Safe, General Purpose, Cotton and our Specialty line featuring the pre-wetted CrushTube™ Swab.

Texwipe Cleanroom Swabs are manufactured with a complete thermal bond, eliminating adhesive contamination\*. Our swabs are manufactured using high-precision automated processes, resulting in swabs constructed to exacting and consistent tolerances, lot coded for traceability and quality control.

### Custom Swabs

At Texwipe® we understand that the best swab for your application may be something unique. We can create a custom swab for your specific applications.

### Sterile Swabs

We can create sterile versions of all Texwipe swabs or something unique and specific to your application.

\*excludes Cotton Swabs


# Cleaning Validation Swabs

Engineered to deliver consistent and accurate sample recovery for cleaning validation for both specific (HPLC, UV-Vis, IMS) and nonspecific (TOC) analytical methods. The swabs are effective with a wide range of diluents.


## Alpha®

Alpha series swab additionally processed to <50 µg/L TOC per swab.


Part No. <b>TX714K</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	4.2mm (0.165")	25.7mm (1.012")	5.2mm (0.205")	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020")	20 swabs per bag 50 bags per case



Part No. <b>TX715</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	4.2mm (0.165")	25.7mm (1.012")	5.2mm (0.205")	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020")	100 swabs per bag 2 inner bags of 50 swabs 10 bags per case




Part No. <b>TX761K</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm (0.268")	2.8mm (0.110")	16.8mm (0.661")	3.2mm (0.126")	3.2mm (0.126")	145.5mm (5.728")	162.3mm (6.390")	20 swabs per bag 50 bags per case



## Absorbond®

Surfactant-free Absorbond Series swab for HPLC/IMS analysis.

Part No. <b>TX716</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	4.2mm (0.165")	25.7mm (1.012")	5.2mm (0.205")	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020")	20 swabs per bag 50 bags per case



# TOC Cleaning Validation Kits

Texwipe's TOC Cleaning Validation Kits are designed to simplify sampling as part of a cleaning validation protocol. The swabs and vials are further cleaned using a proprietary processes to make these kits excellent for use in Total Organic Carbon (TOC) analysis.

The kits are an organized, convenient and secure method for transporting swabs and vials from storage area to sampling sites to laboratory with minimal chance of contamination.

Packaged in cleanroom-compatible, recyclable polypropylene boxes.

Two size options available.



TOC Cleaning Validation Kit

Product Number	Description	Packaging
<b>TX3340</b>	TOC Cleaning Validation Kit - 12 vials + 24 swabs + 12 labels	18 kits/case
<b>TX3342</b>	TOC Cleaning Validation Kit, Bulk - 72 vials + 144 swabs + 72 labels	1 kit/case

# Polyester Swabs

## Microdenier Series

Texwipe's Microdenier Swab Series is made from the cleanest microdenier material, 100% polyester, nylon-free. The complete thermal bond construction eliminates adhesive contamination.

Cleanroom laundered, providing ultra low levels of particles, NVRs (non volatile residues) and ions, and engineered to provide enhanced absorbency and quick wicking performance. Streak-free cleaning with scratch-sensitive surfaces and compatible with QAC (quatarnary ammonium chloride) cleaners and disinfectants. 100% virgin polypropylene handle.


Cleanroom manufactured, made to exacting and consistent tolerances using high-precision automated processes. Lot coded for traceability and quality control.

Packaged in a reclosable, silicone-free bag.



Reclosable, silicone-free bag

Part No. <b>TX714MD</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	4.2mm (0.165")	25.7mm (1.012")	5.2mm (0.205")	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020")	100 swabs per bag 8 bags per case



Part No. <b>TX758MD</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.2mm (0.126")	2.5mm (0.098")	11.0mm (0.433")	2.2mm (0.087")	2.2mm (0.087")	60.0mm (2.362")	70.0mm (2.756")	500 swabs per bag 5 bags per case








Part No. <b>TX761MD</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm (0.268")	3.6mm (0.142")	16.8mm (0.661")	3.2mm (0.126")	3.2mm (0.126")	145.5mm (5.728")	162.5mm (6.398")	100 swabs per bag 8 bags per case



# Polyester Swabs

## Alpha® Series

Texwipe's Alpha® Swab Series is made from the cleanest 100% polyester knit materials with a double layer of fabric for enhanced absorbency. A complete thermal bond construction eliminates adhesive contamination. Cleanroom laundered and manufactured, providing ultra low levels of contaminants. Engineered to resist snagging and abrading. 100% virgin polypropylene handle. Made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.

Part No. <b>TX714A</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	4.2mm (0.165")	25.7mm (1.012")	5.2mm (0.205")	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020")	100 swabs per bag 2 inner bags of 50 swabs 10 bags per case
Part No. <b>TX743B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	4.0mm (0.157")	3.6mm (0.142")	12.0mm (0.472")	3.0mm (0.118")	3.0mm (0.118")	58.3mm (2.295")	70.3mm (2.768")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. <b>TX754B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.5mm (0.138")	3.0mm (0.118")	9.5mm (0.374")	2.4mm (0.095")	2.4mm (0.095")	60.0mm (2.362")	69.5mm (2.736")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. <b>TX758B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.2mm (0.126")	2.5mm (0.098")	11.0mm (0.433")	2.2mm (0.087")	2.2mm (0.087")	60.0mm (2.362")	69.5mm (2.736")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. <b>TX761</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm (0.268")	2.8mm (0.110")	16.8mm (0.661")	3.2mm (0.126")	3.2mm (0.126")	145.5mm (5.728")	162.3mm (6.390")	100 swabs per bag 10 bags per case

## Absorbond® Series

Texwipe's Absorbond® Swabs Series is made from 100% polyester (hydroentangled) nonwoven material double layer for enhanced absorbency. The complete thermal bond construction eliminates adhesive contamination. Cleanroom processed and manufactured, providing ultra low levels of contamination control. Engineered for use on sensitive surfaces. 100% virgin polypropylene handle. Made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.

Part No. <b>TX716</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	4.2mm (0.165")	25.7mm (1.012")	5.2mm (0.205")	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020")	20 swabs per bag 50 bags per case
Part No. <b>TX759B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.2mm (0.126")	1.5mm (0.059")	10.0mm (0.394")	2.2mm (0.087")	2.2mm (0.087")	60.0mm (2.362")	70.0mm (2.756")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. <b>TX762</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm (0.268")	2.6mm (0.102")	16.8mm (0.661")	3.2mm (0.126")	3.2mm (0.126")	145.5mm (5.728")	162.3mm (6.390")	100 swabs per bag 10 bags per case

Sterile available upon request.

# CleanFoam® Swabs

## CleanFoam® Series A

Texwipe's CleanFoam® Swab Series is made with 100% polyurethane foam, 100 pores per inch, for good chemical compatibility with a variety of solutions. The complete thermal bond construction eliminates adhesive contamination.

100% virgin polypropylene handle. Cleanroom processed and made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.

*Foam Series A features 100 swabs per bag.*




Part No. <b>TX706A</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	10.9mm (0.429")	8.5mm (0.335")	27.0mm (1.063")	6.3mm (0.248")	6.3mm (0.248")	54.0mm (2.126")	81.0mm (3.189")	100 swabs per bag 2 inner bags of 50 swabs 10 bags per case



Part No. <b>TX707A</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	15.0mm (0.591")	8.8mm (0.346")	25.0mm (0.984")	6.6mm (0.260")	3.2mm (0.126")	106.0mm (4.173")	131.0mm (5.157")	100 swabs per bag; 2 inner bags of 50 swabs; 10 bags per case



Part No. <b>TX708A</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	19.7mm (0.776")	9.4mm (0.370")	26.4mm (1.039")	6.6mm (0.260")	3.2mm (0.126")	105.8mm (4.165")	132.2mm (5.205")	100 swabs per bag; 2 inner bags of 50 swabs; 10 bags per case



Part No. <b>TX709A</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	4.8mm (0.189")	4.2mm (0.165")	12.0mm (0.472")	4.3mm (0.169")	4.3mm (0.169")	92.5mm (3.642")	104.0mm (4.094")	100 swabs per bag 2 inner bags of 50 swabs 10 bags per case



Part No. <b>TX710A</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	7.0mm (0.276")	5.9mm (0.232")	20.4mm (0.803")	4.3mm (0.169")	4.3mm (0.169")	92.5mm (3.642")	113.0mm (4.449")	100 swabs per bag 2 inner bags of 50 swabs 10 bags per case



Part No. <b>TX712A</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	7.8mm (0.307")	25.7mm (1.012")	5.2mm (0.205")	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020")	100 swabs per bag 2 inner bags of 50 swabs 10 bags per case





# CleanFoam® Swabs

## CleanFoam® Series B

Texwipe's CleanFoam® Swab Series is made with 100% polyurethane foam, 100 pores per inch, for good chemical compatibility with a variety of solutions. The complete thermal bond construction eliminates adhesive contamination.

100% virgin polypropylene handle. Cleanroom processed and made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.

*Foam Series B features 500 swabs per bag.*



Part No. <b>TX740B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.2mm (0.244")	5.0mm (0.197")	17.0mm (0.669")	3.2mm (0.126")	3.2mm (0.126")	146.0mm (5.748")	163.0mm (6.417")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX741B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.4mm (0.134")	3.0mm (0.118")	10.5mm (0.413")	2.4mm (0.094")	2.4mm (0.094")	56.8mm (2.236")	67.3mm (2.650")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX742B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.6mm (0.142")	3.5mm (0.138")	12.0mm (0.472")	3.0mm (0.118")	3.0mm (0.118")	59.0mm (2.323")	71.0mm (2.795")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX751B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.5mm (0.138")	3.5mm (0.138")	16.0mm (0.630")	2.4mm (0.094")	2.4mm (0.094")	51.0mm (2.008")	67.0mm (2.638")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX752B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.4mm (0.134")	3.4mm (0.134")	20.0mm (0.787")	2.5mm (0.098")	2.5mm (0.098")	50.0mm (1.969")	70.0mm (2.756")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX757B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.2mm (0.126")	2.6mm (0.102")	10.0mm (0.395")	2.2mm (0.087")	2.2mm (0.087")	59.7mm (2.350")	69.7mm (2.744")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Sterile available upon request.

# ESD-Safe Swabs


Texwipe's ESD Swab Series features a Stat-Rite® handle, a patented Inherently Dissipative Polymer, to provide clean ESD protection without contributing particulate contaminants to the process or environment. The complete thermal bond construction eliminates adhesive contamination. Packaged in an anti-static bag.

Available in a choice of CleanFoam®, Alpha® Polyester, or Stick.




## CleanFoam® ESD


Part No. <b>TX750E</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.5mm (0.138")	3.5mm (0.138")	10.0mm (0.394")	3.0mm (0.118")	3.0mm (0.118")	78.5mm (3.090")	88.5mm (3.484")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX753E</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.7mm (0.146")	3.4mm (0.134")	10.0mm (0.394")	3.0mm (0.118")	3.0mm (0.118")	80.0mm (3.150")	90.0mm (3.543")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX757E</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.8mm (0.150")	3.4mm (0.134")	10.0mm (0.394")	3.0mm (0.118")	3.0mm (0.118")	80.7mm (3.177")	90.7mm (3.571")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



## Alpha® Polyester ESD

Part No. <b>TX758E</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	4.0mm (0.157")	3.0mm (0.118")	10.0mm (0.394")	3.0mm (0.118")	3.0mm (0.118")	80.5mm (3.169")	90.5mm (3.563")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case



Part No. <b>TX761D</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm (0.268")	4.0mm (0.157")	17.0mm (0.669")	3.0mm (0.118")	3.0mm (0.118")	147.0mm (5.787")	164.0mm (6.457")	100 swabs per bag 10 bags per case



## Stick ESD

Part No. <b>TX769E</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	1.0mm (0.039")	1.0mm (0.039")	9.6mm (0.378")	2.4mm (0.094")	2.4mm (0.094")	56.5mm (2.224")	66.1mm (2.602")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case





# Cotton Swabs

Texwipe's Cotton and Cotton/Foam Swab Series is designed to maximize the characteristics of cotton. Available in 100% cotton or cotton/polyurethane foam. These swabs are bonded securely to a wooden shaft by an aqueous-based adhesive.

Part No. <b>TX700B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	13.5mm (0.531")	13.5mm (0.531")	27.5mm (1.083")	2.5mm (0.098")	2.5mm (0.098")	125.5mm (4.941")	152.5mm (6.004")	500 swabs per bag 10 inner bags of 50 swabs 5 bags per case



Part No. <b>TX705</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	7.0mm (0.276")	7.0mm (0.276")	17.0mm (0.669")	2.5mm (0.098")	2.5mm (0.098")	135.0mm (5.315")	152.0mm (5.984")	1,000 swabs per bag 10 inner bags of 100 swabs 5 bags per case



Part No. <b>TX720B</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	11.8mm (0.465")	10.0mm (0.394")	26.3mm (1.035")	2.5mm (0.098")	2.5mm (0.098")	127.0mm (5.000")	154.0mm (6.063")	500 swabs per bag 10 inner bags of 50 swabs 5 bags per case



# Specialty Swabs & Tools

TX726 is our innovative CrushTube Swab, featuring 91% IPA in an enclosed tube. Packaged in a resealable bag.

Part No. <b>TX726</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	8.3mm (0.327")	n/a	10.0mm (0.394")	9.5mm (0.37")	9.5mm (0.37")	44.0mm (1.73")	47.6mm (1.89")	50 swabs per box 10 boxes per case



TX730 is a 100% nylon swab with micro bristles that grab particles or remove liquids by capillary forces. Ideal for tight spaces. Packaged in a silicone-free, amide-free bag.

Part No. <b>TX730</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	0.8mm (0.031")	0.8mm (0.031")	15.5mm (0.610")	2.8mm (0.110")	2.8mm (0.110")	152.5mm (6.004")	168.0mm (6.614")	50 swabs per bag 10 bags per case



# General Purpose Swabs

Texwipe's General Purpose Swab Series is designed for those customers needing a swab in industrial or general purpose applications. 100% polypropylene handle provides excellent chemical compatibility. Complete thermal bond construction eliminates adhesive contamination.

Packaged in a reclosable, silicone-free bag. Available in Polyester Honeycomb and Foam.



Reclosable, silicone-free bag

## Polyester Honeycomb

Texwipe's General Purpose Polyester Swab Series, made with 100% polyester, was engineered with a honeycomb-type cell designed to maximize absorbency and capture of large particles and dirt in less critical areas.

Part No. <b>TX801</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	13.5mm (0.531")	3.5mm (0.138")	25.2mm (0.992")	5.2mm (0.205")	3.0mm (0.118")	102.3mm (4.028")	127.5mm (5.020")	100 swabs per bag 6 bags per case



Part No. <b>TX802</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.5mm (0.138")	3.0mm (0.118")	9.2mm (0.362")	2.5mm (0.098")	2.5mm (0.098")	60.0mm (2.362")	69.2mm (2.724")	500 swabs per bag 5 bags per case



## Foam

Texwipe's General Purpose Foam Swab Series, made with 100% polyurethane foam, was engineered for applying and absorbing adhesives, greases and liquids in less critical areas.

Part No. <b>TX803</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.6mm (0.142")	3.5mm (0.138")	12.0mm (0.472")	3.0mm (0.118")	3.0mm (0.118")	59.0mm (2.323")	71.0mm (2.795")	500 swabs per bag 3 bags per case



Part No. <b>TX804</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.5mm (0.138")	3.5mm (0.138")	15.5mm (0.610")	2.4mm (0.094")	2.4mm (0.094")	51.0mm (2.008")	66.5mm (2.619")	500 swabs per bag 3 bags per case



Part No. <b>TX805</b>	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	19.7mm (0.776")	9.4mm (0.370")	26.4mm (1.039")	6.6mm (0.260")	3.2mm (0.126")	105.8mm (4.165")	132.2mm (5.205")	80 swabs per bag 5 bags per case



## Swab Solvent Compatibility

## Swab Index

Solvent	Head Substrate		Handle Material	
	Foam	Polyester	Polypropylene	Stat-Rite®
Acetic Acid	U	U		
Acetone	S	U	U	D
Acetonitrile	U	U	U	D
Ammonia, Anhydrous	S	U		
Aniline	S	U		
Benzene	S	U	U	D
Carbon Tetrachloride	S	U		
Chloroform	S	U		
Chlorophenol	D	D		
Chromic Acid	D	D		
Cyclohexane	S	U		
Dimethylformamide	D	S		
Dimethylsulfoxide	D	S		
Ethyl Acetate	S	U		
Ethyl Alcohol	U	U		
Ethyl Chloride	S	U		
Formic Acid	D	U		
Hexane	U	U	U	U
Hydrochloric Acid 48%	D	D	U	U
Hydrogen Peroxide	D	U		
Isopropyl Alcohol	U	U	U	U
Methanol	S	U	U	U
Methyl Ethyl Ketone	S	U		
Methylene Chloride	S	U	U	D
Nitric Acid 100%	D	D	U	D
Nitrobenzene	S	U		
Perchloroethylene	U	U		
Phosphoric Acid 85%	D	U		
Sodium Hydroxide 36%	S	U		
Sulfuric Acid 98%	D	D		
Toluene	S	U	U	D
Trichloroethylene	S	U		
Vinyl Acetate	S	U		
Water	U	U	U	U
Xylene	S	U	U	D

TX700B	Pg. 7
TX705	Pg. 7
TX706A	Pg. 4
TX707A	Pg. 4
TX708A	Pg. 4
TX709A	Pg. 4
TX710A	Pg. 4
TX712A	Pg. 4
TX714A	Pg. 3
TX714K	Pg. 1
TX714MD	Pg. 2
TX715	Pg. 1
TX716	Pg. 3
TX720B	Pg. 7
TX726	Pg. 7
TX730	Pg. 7
TX740B	Pg. 5
TX741B	Pg. 5
TX742B	Pg. 5
TX743B	Pg. 3
TX750E	Pg. 6
TX751B	Pg. 5
TX752B	Pg. 5
TX753E	Pg. 6
TX754B	Pg. 3
TX757B	Pg. 3
TX757E	Pg. 6
TX758B	Pg. 3
TX758E	Pg. 6
TX759B	Pg. 3
TX761K	Pg. 1
TX758MD	Pg. 2
TX761	Pg. 3
TX761D	Pg. 6
TX761MD	Pg. 2
TX762	Pg. 3
TX769E	Pg. 6
TX801	Pg. 8
TX802	Pg. 8
TX803	Pg. 8
TX804	Pg. 8
TX805	Pg. 8

Compatibility is determined at room temperature.

Key: D–Degrades; S–Swells, U–Unaffected



[www.texwipe.com](http://www.texwipe.com)

*For more information please contact us.*

Texwipe Americas	+1 336 996 7046
Texwipe Europe	+31 88 1307 410
Texwipe Asia	+63 049 530 3390
Texwipe China	+65 6468 9433